

## TRANSFORM YOUR MACROSCOPE INTO AN ULTRAMACROSCOPE FOR 3D IMAGING OF CLEARED ORGANS

Maeva Gesson<sup>a</sup>, Mahmoud Kechad<sup>b</sup>, Nicolas Mauclert<sup>c</sup>, Paul Girard<sup>c</sup>, Alain Roussel<sup>c</sup>, Olivier Pierre<sup>d</sup>, Cécile Carvajal<sup>e</sup>, Iris Grosjean<sup>f</sup>, Frédéric Brau<sup>e</sup>

<sup>a</sup>Institut de Biologie de Valrose, Université Nice Sophia-Antipolis, UMR CNRS 7277 - INSERM 1091, 06100 Nice, France

<sup>b</sup>irisPHARMA, Allée Hector Pintus, 06610 La Gaude, France

<sup>c</sup>S2M, Observatoire de la Côte d'Azur, Boulevard de l'Observatoire, CS34229, 06304 Nice Cedex 4, France

<sup>d</sup>Institut Sophia Agrobiotech, Université de Nice-Sophia Antipolis Sophia Antipolis, France.

<sup>e</sup>Institut de Pharmacologie Moléculaire et Cellulaire, Université de Nice Sophia-Antipolis and CNRS UMR CNRS 7275, 660 route des lucioles, 06560 Valbonne, France

<sup>f</sup>Institute for Research on Cancer and Aging of Nice (IRCAN), Université Nice Sophia-Antipolis, UMR7284 CNRS and U1081 INSERM Nice 06107, France

E-mail : brau@ipmc.cnrs.fr

**KEY WORDS** : light sheet microscopy, tissue clearing; whole organ imaging

Serial thin sectioning of whole organs are commonly used to indirectly get the three-dimensional (3D) reconstruction of their histological structures. This approach requires intact sections as well as proper alignment of all the individual images [1]. In the 90's, light-sheet fluorescence imaging combined with clearing methods was firstly used to image whole guinea pig cochlea [2]. More recently, this approach was implemented for bigger samples with a dual illumination path, coupled with the development of clearing reagents [3]. A version of this system is commercially available. Here, we propose an open and versatile light sheet based instrumentation which can easily be implemented on commercial fluorescence microscopes. Many microscopy facilities and labs already possess a detection unit *i.e.* commercial fluorescence microscopes coupled with cameras. Our approach allows a cost effective upgrade of a pre-existing systems into "ultramicroscopes". It consists in a plug and play cylindrical-lens based dual illumination unit with a sample holder mounted on a fast vertical translation stage. Optical fibers coming from laser boxes can directly be plugged on the unit. The two light-sheets can be easily aligned using a specially-designed alignment cuvette and/or with the help of a dynamic overlay of the fluorescence images obtained by a sequential right/left illumination. All the hardware is piloted by the open source MicroManager software [4]. A previous version of the prototype proved its efficiency on glioblastoma imaging in whole mice brains cleared with 3Disco [5]. We now present our last 3D images obtained on this new prototype with 3Disco, iDisco and Ultimate Disco cleared organs.

[1] Ma B, Lin Z, Winkelbach S, Lindenmaier W, Dittmar KE, Automatic registration of serial sections of mouse lymph node by using Image-Reg. *Micron*, Jun; 39 (4):387-96 (2008)

[2] Voie AH, Burns DH, Spelman FA., Orthogonal-plane fluorescence optical sectioning: three-dimensional imaging of macroscopic biological specimens. *J Microsc*, Jun; 1, 70(Pt 3):229-36 (1993)

[3] Dodt HU, Leischner U, Schierloh A, Jährling N, Mauch CP, Deininger K, Deussing JM, Eder M, Ziegglänsberger W, Becker K., Ultramicroscopy: three-dimensional visualization of neuronal networks in the whole mouse brain. *Nat Methods*, Apr; 4 (4):331-6 (2007)

[4] Arthur Edelstein, Nenad Amodaj, Karl Hoover, Ron Vale, and Nico Stuurman, Computer Control of Microscopes Using µManager. *Current Protocols in Molecular Biology* (2010)

[5] Garofalo S, D'Alessandro G, Chece G, Brau F, Maggi L, Rosa A, Porzia A, Mainiero F, Esposito V, Lauro C, Benigni G, Bernardini G, Santoni A, Limatola C. Enriched environment reduces glioma growth through immune and non-immune mechanisms in mice. *Nat Commun*, Mar 30; 6:6623 (2015)